


QuantaVue 4X Real-time PCR



■ Real-Time qPCR System

The quantitative polymerase chain reaction (qPCR) is an established method for the highly sensitive detection and quantification of DNA or RNA. The measuring principle is based on fluorescent signals, which, cycle by cycle, capture the presence of the existing target sequence in real time. The key features of this detection method are outstanding high-performance optics as well as excellent temperature uniformity over 96 samples.

The Unicorn product family guarantees well-founded real-time PCR results as it benefits from peerless temperature control precision in the sample block regardless of the number of samples used. The patented high-performance optics guarantee the outstanding homogeneous excitation and illumination of all individual samples. The Unicorn product family sets new standards of flexibility and precision – for all real-time PCR applications.



■ Redefining Excellence

- Innovative optical detection design - Higher Sensitivity
- Unique real-time scanning technology - Less Optical Cross-talk
- Special hollow-out sample block - Faster & Uniform Thermal Conductivity
- Powerful & Easy-to-Use software - Intuitive Operation, Intelligent Analysis
- Plug & Play - Ex-factory Calibration, Maintenance-free

■ Applications


- Gene Expression Analysis
- Genotyping
- Gene Mutation Analysis
- Pathogen Detection & Quantification
- GMO Detection
- Protein Stability Screening
- miRNA Analysis

Qualified
Precise
Creative
Reliable

Product Features

1. Innovative Optical Design--- Higher Sensitivity, Less Cross-talk


As for the core detector, the Unicorn system further improves the detection sensitivity by using Fresnel lens combined with large area MPPC, a patented new optical detection system. High quality MPPC is used to ensure the detection sensitivity of the system. Meanwhile, it combines the advantages of thin Fresnel lens with short focal length to make the detector closer to the sample, thus effectively reducing the loss of optical signal as well as optical cross-talk between samples, and finally improves the detection sensitivity. Further, the matching of high frequency LED and MPPC makes the detection time of each well and each



The state of the art optics shuttle of the Unicorn system

fluorescent channel very short, so that effectively reducing fluorescence quenching.

In terms of scanning methods, the Unicorn innovatively adopts the time-resolved well-by-well scanning technology. Through the setting of interlaced scanning heads, the problem of cross-talk between wells is significantly reduced from the perspective of spatial distribution. Meanwhile, high-precision scanning heads are used to detect signals of different fluorescent channels in different wells in time sequence (time-resolved), thus eliminating the cross-talk in principle.



Innovative time-resolved well-by-well scanning pattern


- High Sensitivity
- Less Cross-talk
- Fast Scan
- No Edge Effects
- Maintenance-free

2. Outstanding Thermal Cycler --- Superior Uniformity, Precision, Speed

The Unicorn system uses the latest Peltier components through strict screening and testing for excellent performance and stable quality, the unique hollow-out thermal block which has higher ramp rates than standard blocks by reducing the overall mass of metal block, the edge temperature compensation technology based on thermal conductive carbon film and auxiliary heating plate for maintaining tight temperature uniformity throughout the entire sample plate, and the precise temperature control algorithm for fast time to target temperature and faster protocol run times.

As a result, the Unicorn is distinguished by excellent temperature control precision and uniformity of $\pm 0.2\text{ }^{\circ}\text{C}$ (over 96 wells) and first-class heating rates of up to $6\text{ }^{\circ}\text{C/s}$. In addition, to ensure the highest specificity for different assays, the device is equipped with a gradient function with a spread

of up to $36\text{ }^{\circ}\text{C}$. This combination makes the entire system absolutely ideal for any real-time PCR application.




Special hollow-out thermal block for reliable, consistent results

- Compatible with 0.1ml or 0.2ml Low Profile 96-well plate, 8-strip tube, single tube(transparent, frosted, and milky white are applicable)
- Precise temperature control of $\pm 0.2\text{ }^{\circ}\text{C}$ across the entire sample block
- Innovative thermal block with outstanding ramping rates of up to $6\text{ }^{\circ}\text{C/s}$
- Gradient function over 12 columns with a $36\text{ }^{\circ}\text{C}$ spread

3. Impressive performance you can trust to generate high-quality data



The combination of innovative optics and unique thermal cycler with incomparably precise temperature control ensures ideal amplification results. This makes the Unicorn a trustworthy partner for quantitative real-time PCR applications.

- Excellent reproducibility and 10-log dynamic range
- Precise quantification with 1.33-fold discrimination
- True Five-Target Multiplexing




Excellent reproducibility for reliable results.

Amplification of a plasmid template in 96 wells. The mean ct value of 11.86 with a standard deviation of 0.05 was determined automatically, illustrating highly homogeneous amplification results obtained by the Unicorn.





Broad linear dynamic range ensuring accurate quantification.

The example of plasmid DNA amplification shows an optimal linearity over 10 orders of magnitude from 2.5 ng/μl to 0.12 fg/μl, illustrating the broad linear dynamic range of the system




Sensitive detection and high-confidence target discrimination down to 1.33-fold.

The results of amplification with different amounts of plasmid were obtained by Taqman method (the template concentration was 500 Pg/μl, 250Pg/μl, 125Pg/μl, 30Pg/μl, 20Pg/μl, 13.3Pg/μl in turn from left to right). Graph shows 4 replicates for each dilution. The local enlargement map shows that 2 fold, 1.5 fold and 1.33 fold differences were clearly distinguished.



Superior Uniformity, Precision and Speed.

The Unicorn thermal cycler exhibits high average ramp rates, rapid setting time, and tight thermal uniformity. This graph shows the temperature measured by probes in 20 wells across a sample block. The traces are nearly indistinguishable due to the tight uniformity with $<\pm 0.2^{\circ}\text{C}$ temperature fluctuation.



Multiplexing for richer datasets.

The Unicorn system can discriminate up to five targets in a single reaction well. The fluorescent data from four FAM-HEX-ROX-Cy5 channels of four different target genes in a single well showing there is almost no cross-talk between different channels.

4. Special Analyzer Software

Unicorn Analyzer Software is powerful software for a powerful instrument. It accommodates user needs and different types of experiments with intuitive navigation and customizable settings. The logical, clear arrangement of all the tools, the intuitive handling, and last but not least the parameter-oriented memory and programming concept make the software easier to use.

- User-friendly and clearly structured
- Intuitive navigation and customizable settings
- Integrated analysis algorithms with automatic analysis
- Comprehensive analysis modules for multiple applications
- License-free, free updates



Intuitive navigation and customizable settings

To further speed up your process, get the experiment started with the quick plate feature, and then enter or edit well information on your own schedule - before, during, or after the run has finished.




Intelligent Analysis

To ensure that operation is as simple as possible, many of the steps occur automatically – such as performing normalized gene expression analysis (ddCt), you take your results a step further without having to export data and set up your own analysis macro.





Personalized Data Export

The software allows you to save predefined analysis settings for auto-exporting run data into their format of choice, including Excel, PDF, txt export format.




Absolute quantification

The module enables analysis of genes of interest with the use of a standard curve. Additional flexibility is achieved by importing standard curves from other experiments.




Relative quantification

With this module, you can customize groupings of data within projects for a thorough comparison of data. The module also includes integrated correlation, volcano, and cluster plot analysis, with the ability to drill down to amplification plots..




5. Comprehensive Function

The software has built-in data analysis modules with automatic baseline subtraction and threshold calculation for determining Ct values or possible standard curves and PCR efficiencies. Further analyses can likewise be conducted automatically, such as absolute or relative quantifications. The software also includes analysis methods for probe-based allelic discrimination and the use of a positive/negative analysis via the end-point detection of samples.




Positive/Negative Analysis.

The figure shows that the positive and negative African swine fever virus samples are well differentiated.



Melt curve analysis.

In this experiment, 96 replicates of human genomic DNA were amplified followed by a dissociation step with a melting temperature (T_m) of 77.9°C (SD 0.07°C), showing good specificity of the reaction.



Protein Stability Screening.

The figure shows that the protein exhibits different T_m values when binding with ligands of different concentrations (1mM, 0.1mM, 0mM), which indicates that the protein thermal stability changed.

■ Specification

Model	QuantaVue 4X	
Block capacity	96	
Sample volume	10-50µl	
Heating/cooling method	Peltier	
Maximum ramp rate	6 °C/s	
Temperature setting range	4-100 °C	
Heated lid	Electronic automatic lid	
Temperature accuracy	±0.2 °C	
Temperature uniformity	±0.2 °C	
Gradient zone	12 columns	
Gradient range	1-36 °C	
Excitation source	Long-life, high performance LEDs	
Detector	Highly sensitive MPPC with Fresnel lens	
Scanning principle	Time-resolved scanning technology	
Detector position	Top of the block	
Excitation/detection range	455-650nm/510-715nm	
Fluorescence channel	4 channels	
Detection sensitivity	1 copy of the target sequence	
System sensitivity	Detect differences as small as 1.33-fold in target quantities in singleplex reactions	
Dynamic range	10 orders of magnitude	
Dye compatibility	FAM/SYBR Green, VIC/JOE/HEX/TET, NED/TAMRA/Cy3, JUN, ROX/Texas Red, Mustang Purple, Cy5/LIZ	
Data Analysis Modes		
Absolute quantification	Relative quantification	Endpoint qualitative analysis
Melt curve analysis	Protein Stability Screening	Genotyping
Data Export		
Customizable reports containing run setting, data graphs, and spreadsheets can be directly exported or saved as Excel, txt, PDFs		

Shanghai Unicorn Medical Technology Co., Ltd.

Room J, Building 4, B Zone, No. 925 Yecheng Road, Jiading Industrial Area, Shanghai

Tel: +86 (021) 6082 8473 Cel: +86 173 1710 8636 Email: info@unicornlifescience.com

www.unicornlifescience.com

Information is subject to change and/or updating without notice
@ 2023 Unicorn. LIFE-REALTIMEPCR-202305