




## Innovative Temperature Control

Maintain a constant temperature of 22°C

### Temp. Control Technology


The unit features a heating device at the bottom and a refrigeration evaporator at the top. During the cooling process, the cold air sinks, while during the heating process, the hot air rises. This design ensures rapid heating or cooling, regardless of the process involved.



Air Circulation

### Ultimate Uniformity

The cabinet maintains temperatures within  $\pm 1^\circ\text{C}$  throughout the entire compartment, ensuring confidence that items are stored at the correct temperature regardless of their placement in the chamber.



### Optimized Fluctuation

Maintain superior temperature stability @  $22\pm 2^\circ\text{C}$ , avoiding rapid and significant changes in temperature. This ensures that the stored platelets are consistently stored in their optimal environment.

### Fast Recovery

Faster temperature recovery even after prolonged door openings. The forced-air circulation system maintains consistent temperatures, ensuring that contents are stored at the right temperature even when the unit is opened multiple times per hour.

### Precision Controller

- A digital microprocessor temperature controller with LCD display and alarm back up power.
- Accurately controlling the cooling or heating process to maintain a constant temperature of  $22^\circ\text{C}$ .
- Temperature setting range:  $20^\circ\text{C} \sim 24^\circ\text{C}$ .
- The precision of display and control is  $0.1^\circ\text{C}$ .
- Alarms: High/low temperature, power failure, sensor error, door ajar, backup battery for 72 hours.



C309 controller



Heating status



Cooling status



# Platelet Incubator

Platelet storage @ 22°C ± 2°C



SDC-5



SDC-10

## Ultimate insulation

The cabinet is insulated with a thick layer of high-density PU foam, effectively sealing and maintaining the internal temperature. The door, constructed with the same insulation material, is equipped with double-layer vacuum glass coated with LOW-E film that ensures secure platelet preservation while providing enhanced safety measures.



Insulated foamed glass door

## Electrically heated glass door


Fog or condensation often occurs when opening and closing the door. However, this can be easily prevented with the electrically heated glass door, which quickly evaporates any moisture.



Complete SS interior

## High quality materials

The units feature a high-quality coated cold-rolled steel exterior and stainless steel interior, offering antibacterial, anti-corrosive, and easy-to-clean properties. The agitator and racks are also constructed with stainless steel, ensuring reliability and hygiene.



SDC-10

## System composition

The system consists of an agitator, stainless steel racks, and stainless steel drawers. The stainless steel racks are securely attached to the agitator platform and move in sync with the agitation motion. This design helps maintain platelets in a suspended state, preventing clumping or aggregation.

## Ultra-long lifespan motors

Imported coaxial motors offer numerous advantages such as low heat generation, low power consumption, minimal noise levels, and a high reliability.

## Gentle glide rails

Imported SUS 304 stainless steel ball bearing sliders have a small gap, low noise, and can continuously operate with a maximum load of 50 kg.





# Platelet Incubator

## Specifications



	Model	SDC-5	SDC-7	SDC-10	SDC-15
Specifications	Drawer QTY	5	7	10	15
	Cooling system	No Frost	No Frost	No Frost	No Frost
	Defrost	Auto	Auto	Auto	Auto
	Refrigerant	R134a	R134a	R134a	R134a
	Cooling system	Air Cooling	Air Cooling	Air Cooling	Air Cooling
	Noise(db)	53	53	53	53
	Ambient temperature (°C)	10~32°C	10~32°C	10~32°C	10~32°C
	Temperature range (°C)	22°C ±2°C	22°C ±2°C	22°C ±2°C	22°C ±2°C
	Compressor QTY	1	1	1	1
	Sensor	NTC	NTC	NTC	NTC
Cooling system	Temperature controller	Microprocessor	Microprocessor	Microprocessor	Microprocessor
	Display	LCD display	LCD display	LCD display	LCD display
	Voltage(V)	187~242	187~242	187~242	187~242
	Frequency(Hz)	50Hz	50Hz	50Hz	50Hz
	Power Consumption(kWh/24h)	2.7	2.7	4.4	4.4
	Power (W)	197W	197W	280W	280W
	(Current(A)	1.3	1.3	1.5	1.5
	External material	Coated steel	Coated steel	Coated steel	Coated steel
	Internal material	SUS 304	SUS 304	SUS 304	SUS 304
	Insulation	CFC-free PURF	CFC-free PURF	CFC-free PURF	CFC-free PURF
Power	Storage Capacity (450ml blood bags)	10	14	20	30
	NT/GT(kg)	65/87	68/90	85/109	92/115
	Exterior size(W*D*H)(mm)	522*600*1056	522*600*1056	522*600*1305	522*600*1305
	Interior size(W*D*H)(mm)	427*360*467	427*360*467	432*462*728	432*462*728
	Package size(W*D*H)(mm)	590*710*1200	590*710*1200	590*710*1450	590*710*1450
	Shipping CBM	0.5	0.5	0.61	0.61
	20GP/40GP/40HQ	24/48/96	24/48/96	27/60/60	27/60/60
	High/low temperature	Y	Y	Y	Y
	Power failure	Y	Y	Y	Y
	Controller failure	Y	Y	Y	Y
Materials	Sensor error	Y	Y	Y	Y
	Door ajar	Y	Y	Y	Y
	Power failure alarm (For controller)	72 hours	72 hours	72 hours	72 hours
	Castor	Y	Y	Y	Y
	levelling feet	N/A	N/A	N/A	N/A
	Access port	Optional	Optional	Optional	Optional
	USB Interface	N/A	N/A	N/A	N/A
	Illumination	UV lamp	UV lamp	UV lamp	UV lamp
	Thermal printers	Y	Y	Y	Y
	Remote alarm port	N/A	N/A	N/A	N/A
Measurement	Door type/Qty	Foamed glass/1	Foamed glass/1	Foamed glass/1	Foamed glass/1
	Data logger	Optional	Optional	Optional	Optional
Alarms					
Accessories					



	Model	SDC-5	SDC-7	SDC-10	SDC-15
Capacity	Rack layer	5	7	10	15
	Rack material	SUS 304	SUS 304	SUS 304	SUS 304
	Drawer QTY	5	7	10	15
	Drawer material	SUS 304	SUS 304	SUS 304	SUS 304
	Platelet bag(450ml) Qty	10	14	20	30